



Shown CK8C50 model

## Features:

- Simple, user-friendly operation
- Operates in 10 seconds or less
- Guaranteed accuracy is 2% within the range of 50% to 100% nominal volume
- Can use a range of test solutions with dye pellets (available for purchase)
- Works with all liquid handler robots that use well plates, multi-channel pipettes, and electronic pipettes
- Compatible design for seamless programming in liquid handlers
- Easy to read markings on a white polystyrene cartridge
- Measures up to 8 channels simultaneously
- Dye is included with the cartridge (dye is non-toxic and acceptable for drain-disposal)
- Glass capillary tubes minimize evaporation of dye prior to reading results
- Enables validations between calibrations and for testing and calibrating with liquids representative of real samples



## Graduations:

Item #	Product Name	Minimum	Maximum	Major Tick	Minor Tick	Working Range
CK8A1-6	Checkit Go 1 $\mu$ L	0.4 $\mu$ L	1.2 $\mu$ L	0.1 $\mu$ L	0.02 $\mu$ L	0.5 – 1 $\mu$ L
CK802-6	Checkit Go 2 $\mu$ L	0.8 $\mu$ L	2.5 $\mu$ L	0.5 $\mu$ L	0.05 $\mu$ L	1 – 2 $\mu$ L
CK8V5-6	Checkit Go 5 $\mu$ L	2 $\mu$ L	6 $\mu$ L	0.5 $\mu$ L	0.1 $\mu$ L	2.5 – 5 $\mu$ L
CK8J10-6	Checkit Go 10 $\mu$ L	4 $\mu$ L	12 $\mu$ L	1 $\mu$ L	0.2 $\mu$ L	5 – 10 $\mu$ L
CK8K20-6	Checkit Go 20 $\mu$ L	8 $\mu$ L	22 $\mu$ L	1 $\mu$ L	0.5 $\mu$ L	10 – 20 $\mu$ L
CK8C50-6	Checkit Go 50 $\mu$ L	20 $\mu$ L	55 $\mu$ L	5 $\mu$ L	1 $\mu$ L	25 – 50 $\mu$ L

## Shelf life:

Item #	Product Name	Shelf life in months
CK8A1-6	Checkit Go 1uL	16 months
CK802-6	Checkit Go 2uL	16 months
CK8V5-6	Checkit Go 5uL	16 months
CK8J10-6	Checkit Go 10uL	16 months
CK8K20-6	Checkit Go 20uL	12 months
CK8C50-6	Checkit Go 50uL	9 months

## Checkit App:

Checkit App creates thorough test reports which can be viewed on your phone or shared to your preferred location.



Next Advance, Inc., 2113 NY 7, Troy, NY 12180, USA  
Phone: +1 (518) 674-3510 ; [info@nextadvance.com](mailto:info@nextadvance.com)  
[www.nextadvance.com](http://www.nextadvance.com)